

REPRO RESULTS gives you options!

Optimize reproductive efficiency with the IDEXX Milk Pregnancy Test.

IDENTIFY OPEN COWS ACCURATELY & EASILY

- As early as 28 days postbreeding and day 60 postcalving
- Expands the use of existing milk quality samples
- Accurate confirmation of pregnancy status
- Non-invasive means of open diagnosis
- Won't disrupt cow management
- Reduces lock-up time
- Shortens regular vet exams
- Ideal for montly reconfirms & dry checks



CONFIRM PREGNANCY WITH CONFIDENCE

- Lab based test uses trusted ELISA technology to detect pregnancy-associated glycoproteins (PAGs)
- Provides accuracy on par with ultrasound/palpation
- 99% Accuracy diagnosing open cows
- Clear and definitive results with same day report turn around
- Upload data directly to Dairy Comp

| Pregnant cows | | IDEXX Milk Test | | | | Performance | |
|---------------------|-------|-----------------|------|---------|------|--------------|---------|
| | | Total | Preg | Recheck | Open | Sensitivity* | Recheck |
| Days After Breeding | 35-75 | 381 | 331 | 42 | 8 | 97.6% | 11% |
| | >75 | 542 | 538 | 1 | 3 | 99.4% | 0.2% |

*Sensitivity excludes IDEXX Milk recheck results.